

CHEM 8M Workbook Organic Chemistry Lab II







	\sim
	Í
\sim	Į į

OH

CONTENTS		
١.	SAFETY RULES	
II.	Writing Guidelines: Experimental Methods and Sample from Caitlin's Thesis	
III.	Safety Scavenger Hunt & Lab Basics Worksheets	
IV.	** Lecture Note Templates **	
V.	NMR Problem Sets	
VI.	IR and NMR Reference Tables	
VII.	Functional Group Table	

Caitlin Binder, PhD University of California, Santa Cruz Department of Chemistry and Biochemistry

Welcome to the CHEM 8M Workbook!

The 8M Workbook is designed to set you up for success in the organic chemistry lab and help you understand our experiments this quarter. Use this workbook for note-taking during *weekly lab lectures* and for *foundational lab activities* such as safety and writing. *Do you find it challenging to take notes while listening to an instructor*? The *lecture note templates* in this workbook are intended to let you process the material in class instead of just writing letters and symbols to 'figure out later!'

This **workbook and its contents are available for free on Canvas** for you to print out on your own, copy by hand, or download to a tablet before class. Purchase of the **8M Workbook** and/or **Lab Manual** from the UCSC Copy Center is optional. The benefit is to have all the paper 3-hole punched copies that are 'binder-ready' (pun unavoidable :P). This workbook could also be used re-write notes (learning hack!).



I look forward to passing on my organic chemistry knowledge and learning from you!

Sincerely, Caitlin Binder, Ph.D. Follow me on Instagram! @AcroChemist

Title	& Order	Description When to use it	
VIII.	SAFETY RULES	Read and bring to the first lab meeting.	
IX.	Safety Scavenger Hunt & Lab Basics Worksheets	Embark on a 'treasure hunt' to find lab safety items in your new workspace © Then, revisit fundamental techniques with common lab equipment.	bring to the first lab meeting.
Χ.	Writing Guidelines	general style and tips for technical writing; guidelines for abstracts	Read for first lab report; bring to lab, Exp 3
XI.	** Lecture Note Templates	Titles, outlines, structures, and example problems are given to help you stay on track and make it fun!	** bring to Friday lab lecture each week **
XII.	NMR Problem Sets	Practice problems on nuclear magnetic resonance (NMR) spectroscopy	bring to lab on the day of Exp 0
XIII.	NMR Tables of Values	Reference for predicting NMR spectra for organic molecules	bring to lab for Exp 0 and each followed lab/lecture

CONTENTS

The 8M Workbook does NOT contain...

- The course syllabus, including
 - o lab & lecture schedule with assignment due dates see Canvas
- Lab Manual lab PDFs with background, procedures, pre- and post-lab questions for each lab are
 posted on Canvas with the option to order a printed Lab Manual from Bookstore see Canvas

Other required 8M materials...

- Lab Notebook any bound, designated notebook for your experiments; duplicate pages NOT required
- **Calculator** bring to lab each day; no cell phones used in the lab as calculators

Provided, and shared materials in the lab...

- Safety goggles
- Lab Coat and disposable nitrile Gloves

LABORATORY SAFETY RULES

Safety First!

The instructors wish students an enjoyable lab experience! The best way to set you up for success is to provide you with lab common sense and etiquette, as much of the surroundings and techniques will be new to you. In general, be respectful of the space and follow instructions and you'll be fine, however, we may need to enforce repeated violations of the rules with point penalties, either for an individual or the whole section, but we hope it won't come to that! Thanks for reading and doing your best as you learn.

1. **Safety goggles must be worn** at all times when anyone in the room is working with chemicals, especially yourself! You are welcome to step outside to defog goggles if necessary.

2. NO food, drinks, or gum are allowed anywhere in the labs or in your mouth while you're in the labs. All water bottles, snacks, etc. should be left on the table outside the lab, not hidden in your bag. You may step into the hallway for a guick drink or snack during appropriate times, as long as you're not neglecting the experiment or your partner. Please check with your TA beforehand.

- 3. Appropriate lab attire must be worn at every lab. Students cannot go home to change.
 - OK LAB ATTIRE: Pants or long skirt, short or long-sleeve shirt, closed-toe shoes that cover the entire top of the foot. Long hair and loose clothing are confined or tied back.
 - NOT OK: Shorts or short skirts (no exposed ankles), leggings/tights, cropped pants that expose ankles, ripped pants that expose skin, tank tops, sandals, ballet flats, or any other shoes that expose the tops of the feet (Crocs and Tom's are NOT OK!). High heels, baggy clothing, and dangling jewelry are strongly discouraged.

4. Lab coats must be worn over appropriate lab attire (see above).

5. NO running, fighting, or other acts of mischief.

6. **NO visitors**, including pets and side-kicks.

7. Know the **locations of emergency equipment** including fire alarms, fire extinguishers, chemical fume hoods, safety showers, and emergency eye washes.

8. Notify your instructor immediately of any injury, spill, fire, or explosion. You may clean up small spills (less than a few milliliters) yourself, but let the TA know. You're not in trouble unless you do it on purpose!

9. Keep your lab space clean and organized throughout the experiment. Backpacks, purses, jackets, phones, etc. are not allowed where chemicals are being used.

10. Never leave an ongoing experiment unattended. If you need to leave the room, be sure a neighbor is watching your experiment.

11. Unless otherwise specified, dispose of broken glassware in broken glassware boxes only, including ceramics and disposable glass pipets. NO paper or other items in the broken glass boxes. NO PIPETS OR OTHER GLASSWARE IN THE TRASH! That's not cool to the staff and you'll lose points.

12. DO NOT TASTE ANYTHING IN THE LAB. EVER.

13. Never remove chemicals or equipment from the labs or stockroom without permission.

14. **NO unauthorized experiments**. Stick to the given procedure.

15. Follow appropriate procedures for inserting glass into a stopper and/or have the stockroom or your TA assist you. Seriously, students have hurt themselves by not paying attention.

16. Wash your hands and arms with soap and water before you leave the lab, no matter what!

17. Always know the **hazards** as well as the physical and chemical properties of the materials used. Your lab notebook should include a brief note on the safety hazards for each chemical being used based on the safety table within in the experiment.

18. **Read labels carefully and know the name of chemicals with which you are working**. Read labels twice. Read labels twice.

19. Label all containers with chemical/mixture names, your name, and the date before anything goes into that container.

20. Use pluringes and pipet bulbs with glass pipets. NEVER pipet by mouth. It's gross.

21. Check all **glassware for cracks and cleanliness** before using...or you'll be sorry later that you didn't.

22. **Avoid contamination**. Take only what you need from reagent bottles and NEVER return unused chemicals to the original bottle that other students are sharing.

23. **Fume hoods** are often used to minimize chemical exposure. Handle chemicals six inches into the hood. DO NOT PUT YOUR HEAD IN THE HOOD and DO NOT KNEEL IN FRONT OF THE HOOD, or anywhere in the lab.

24. Wash all glassware before leaving lab for the day and return all shared equipment to the designated space. DO NOT WEAR GLOVES WHILE WASHING GLASSWARE.

25. **Dispose of all waste as instructed in the lab handout or by the TA.** Read waste container labels carefully to be sure it's going to the right place. Waste containers are typically in the fume hoods. Let your TA know if a waste container is full. DO NOT LET THE WASTE CONTAINERS OVERFLOW! *Seriously, who does that*?!

26. NO use of flame in the lab. Nearly everything in the organic chemistry labs is flammable.

27. *Wear gloves* when appropriate in the lab and *change your gloves* if you get chemicals on them. <u>They're cheap</u>! Gloves are only a first line of protection. They do not make you invincible! Take off gloves and wash your hands before you leave the room. **DO NOT touch door handles, your face, computers, or phones with gloved hands.**

28. *Minimize chemical exposure* (ex. keep containers capped) and treat every chemical as if it were hazardous.

29. No cell phones or electronic devices are allowed to be used in the labs. If you'd like to take a picture or video of your experiment, ask your TA for permission, but take your gloves off first.

30. **Treat your TA and labmates respectfully.** Please adhere to your TA's instructions and additional safety guidelines announced. Reach out to instructors if you are having interpersonal issues in the lab and we'll help as much as we can.

31. Please restore your locker to its original status. All equipment must be clean and organized in the drawer, whether or not you used those items.

- Check the equipment list and reference the 'perfect drawer' pic on the bulletin board in the lab.
- Replace any missing or broken items with a quick trip the stockroom. Do not bring broken items with you place them in the glass or solid waste instead.
- Check for missing, dirty, broken, and extra items that may have fallen into your drawer.

Writing Guidelines: Experimentals & References

Please review **Parts A-B** of the writing guidelines from CHEM 8L (available on 8M Canvas site) for general format and writing style (passive voice, no personal pronouns).

Part D: Experimental Methods and Compound Characterization

Experimental methods and compound characterization are found at the end of scientific journal articles, dissertations, and other technical documents. This gives the reader instructions on how to recreate the experiment and confirm the structure of the newly synthesized compounds. The format and general content differs depending on the field. Students will include this section at the end of several lab reports using the generally accepted guidelines followed by synthetic organic chemists: one <u>General Methods</u> paragraph followed by <u>one additional paragraph per compound</u> synthesized. *A sample Experimental Methods section is attached and contains much more information than CHEM 8M students are expected to include.* Use passive voice and past tense.

General Methods

Reagents and by-products do not get full descriptions but are mentioned in the "General Methods" section with the following statement: "**All reagents were commercially available**, unless otherwise stated." Typically researchers would then describe how reagents and solvents were purified, but *this does not apply to 8M students*. Next, list the **specifications for IR** (medium for analysis, such as salt plates or Teflon) only if used in the experiment.

Experimental Methods & Characterization

Following general methods, each organic compound or reaction gets its own paragraph (one paragraph per reaction/compound). Some or all of the following should be included in the experimental methods and compound characterization section. This is based on experimental techniques students utilized in the lab.

• Reaction scheme - including reactants, reagents, products, solvent(s), and % yield

(structures and reaction schemes can be hand-written)

- Full chemical name of product in bold (common and/or IUPAC)
- Brief description of reaction set up and workup including...
 - Names and amounts of each reactant and reagent (mmol and mL or mg)
 - Name and amount of solvent (mL)
 - o Order of addition, if pertinent, and reaction conditions (time, temperature)
 - o Description, name, and amount of product obtained and % yield:
 - Ex. "Benzhydrol was obtained as a clear liquid (1.00 g, 87% yield)."

Characterization follows in the same paragraph (after reporting the yield) and includes some or all of the following. There is no characterization for the dye lab.

- Melting point or boiling point
- Distinctive IR stretch(es) one or two distinguishing peaks, such as carbonyl or O-H stretches

Sample from Binder, C. M. "Novel (-)-β-Pinene-Derived Amino Alcohols as Asymmetric Directors for the Addition of Organozinc Reagents to Aldehydes" *UC Santa Cruz*, **2010**.

EXPERIMENTAL METHODS AND COMPOUND CHARACTERIZATION

General Methods.

All reagents were commercially available, unless otherwise stated. All air and moisture sensitive reactions were carried out under argon atmosphere using flame- or oven-dried glassware and standard syringe technique. Tetrahydrofuran (THF), dichloromethane (DCM), cyclohexane, triethylamine (Et₃N), morpholine, *tert*-butanol (*t*-BuOH), and dimethyl sulfoxide (DMSO) were distilled over CaH₂. Oxalyl chloride was distilled without drying agent prior to use. Column chromatography was carried out with Silica Gel 60. IR spectra were carried out on NaCl plates with v_{max} in inverse centimeters. Optical rotations were obtained on a digital polarimeter at 20 °C. High resolution mass measurements were obtained on a benchtop ESITOF mass spectrometer.



(+)-Nopinone. NalO₄ (44.96 g, 210 mmol) was added to a 2-L round-bottom flask equipped with a magnetic stir bar and dissolved in water (300 mL), CCl₄ (200 mL), and CH₃CN (200 mL). (–)- β -Pinene (13.88 g, 102.0 mmol)

was added followed by RuCl₃-3H₂O (457 mg, 1.7 mmol). The reaction was stirred overnight while open to the atmosphere (24 h). The crude reaction mixture was filtered through a pad of celite and rinsed with DCM, creating two distinct layers. The aqueous layer was extracted with DCM (3 x 100 mL). The combined organic extracts were washed with water (2 x 30 mL), dried (MgSO₄), filtered, and concentrated *in vacuo* to a black liquid. This was purified by column chromatography (500 mL SiO₂, 100% hexane to elute β -pinene, 4:1 Hexane/EtOAc to elute nopinone) and the nopinone fractions were concentrated to a clear oil (8.3 g, 59% yield). bp 74-76 °C (2 mm Hg), [α]_D²² +34.43° (*c* 4, MeOH), IR (neat) 1714 cm⁻¹.

Part E. Format for Literature References (for future reference in lab reports)

There is a standard A.C.S. (American Chemical Society) format for listing references in the chemical literature that you are required to follow (<u>http://pubs.acs.org/books/references.shtml</u>). This format, illustrated below, must be used in the reference section of your report, if appropriate. Be sure to document all assertions and past work described in your reports with a footnote. Footnotes can be referred to more than once. Use superscripts with corresponding numbered references at the bottom of the page or at the end of the report.

Details matter!

- Pay special attention to punctuation (period / comma / semicolon).
- Years in bold,
- Volumes and titles in italics.

BOOKS

Author's last name, first initial, *Title of Book,* Publisher: City of publication, **Year of pub**.; pages used.

Examples

Crews, P.; Rodríguez, J.; Jaspars, M. Organic Structure Analysis, 2nd Ed.; Oxford: New York, **2010**; pp. 67-70.

Palleros, D.R., *Experimental Organic Chemistry;* Wiley: New York, **2000**; pp. 61-70.

JOURNALS

Author's last name, initials.; 2nd author's last name, initials.; (continue for each author). *Journal abbrev.* **Year**, *Vol.,* first to last page of article.

*Proper journal abbreviation used in italics, **year in bold**, volume in italics, no issue number

Examples

Tansakul, C.; Lilie, E.; Walter, E. D.; Rivera III, F.; Wolcott, A.; Zhang, J. Z.; Millhauser, G. L.; R. Braslau, R. *J. Phys. Chem. C*, **2010**, *114*, 7793-7805.

Sanchez, L. M.; Lopez, D.; Vesely, B. A.; Della Togna, G.; Gerwick, W. H.; Kyle, D. E.; Linington, R. G. J. Med. Chem., **2010**, 53, 4187-97.

Woehrmann, M. H., Gassner, N. C., Bray, W. M.; Stuart, J. M.; Lokey, S. J. Biomol. Screen. 2010, 15, 196-205.

WEB SITES

Use full website addresses to allow the reader to locate referenced material on the web. **Be wary of the content. The info on the web is usually not peer reviewed, and can be erroneous!** If you do cite a website, include the date the website was accessed.

Example

http://organicchemistry.wordpress.com/2007/08/18/tips-for-writing-organic-chemistry-lab-reports/ accessed 7-23-09.

Safety Scavenger Hunt

Welcome to the organic teaching labs! We want your time in the lab to be just as enjoyable as it is educational. You'll learn how to handle equipment properly to eliminate or minimize contamination and chemical exposure. Please keep your lab space neat and organized at all times, similar to cleaning the kitchen dishes as you're cooking or baking. Instrumentation is shared with the rest of the class so we'll work together to stay safe.

Student lockers, reagents, equipment, and instrumentation are organized, set up, and otherwise maintained by the stockroom staff, all of which are very nice people. Please value their time, efforts, and kindness by **leaving the lab as you found it or better**.

Lab Orientation Content

- A. Lab Safety Scavenger Hunt
- B. Equipment Locker

Part A. Lab Safety Scavenger Hunt

Get to know your space!

- Work with your lab-mates to find the following items in the lab and their corresponding tag.
- Make a lab map on the next page and mark the locations by number.
- Fill in the blanks in this packet with the information on the tags.
- There are some items for which there are multiple locations, such as sinks, but only one tag.
- Other items may not be in the room at all!

Er	mergency Response	<u>Day-to-Day</u>	Equipment
1.	Fire Extinguisher	9. Balance Station	16. Equipment Room (GC & IR)
2.	Fire Alarm	10. Sink	17 . Rota-vap
3.	Safety Shower	11. Chemical Waste Station	18. Water Re-circulation Pumps
4.	Eyewash Station	12. Dry Waste Box	(water lines)
5.	Evacuation Procedure	13. Chemical Fume Hoods	19. Ring stands
6.	First Aid Kit	14. Reagent Station (Chemical	20 . Clamps
7.	Broken Glassware Box, Dust	<u>Re</u> acting <u>Agents</u>)	21. Vacuum Tubing
	Pan & Broom	15. Disposable Gloves	22. Hot/stir plates
8	Spill Control Center		

Miscellaneous

23. Your TA – go say hi!

24. One-word hazard definitions & precautions

25. <u>NFPA Labels</u> - Copy and color the NFPA label description from the bulletin board then classify the sample labels posted.

26. Lab coats

LAB MAP, Thimann Labs

Get familiar with the new lab space by adding the locations of items to make a map by number (1-26).

Keep this map handy for future reference!

Start by adding the DOORS and WINDOWS for reference.

HALLWAY



- 1. Fire Alarm: find the closest one in the hallway
- * Located at the _____building entrance and exit left of the elevator.
- * Pull this alarm only if you see a fire. Don't assume someone else has called it in.
- * Alert by-standers by yelling "_____!"
- * Individuals should also notify the fire department by calling ______.

<u>2. Fire Extinguisher</u>: find the closest one in the hallway

* Located in the hallway between rooms _____ and _____.

* Report the fire, _____, and EXIT the building.

* Only ______ should attempt to extinguish a fire

(...so you should probably not be using this, but it's good to know where it is).

3. Safety Shower

* To be used if student is splashed with a considerable amount of ______that cannot be

* DO NOT pull it to test it! Only do it if you need it. It puts out a set, large amount of water.

* If needed, disrobe and stand under the shower for ______ to wash away the chemical and reduce contact (______).

* Call 911 for severe cases.

4. Eyewash Station

* You shouldn't need this because you should be _____!

* Hold the eyelid open in the running eyebath for _____ min.

* Call ______ or go to the ______ if the injury requires further medical attention.

5. Evacuation Procedure

* Take the ______ to evacuate the building instead of the ______.

* Follow the map to the rendezvous point and take a quick spin on the _____.

* Follow the path down the little hill, cross the street, and take a selfie next to the 'Thimann Laboratories

Emergency Assembly Area' sign. Show your TA for credit!

6. First Aid Kit

* Located in the	* Sufficient for	&		·	
* All injuries must be reported to the TA followed by the completion of an "					
Form" attained from the stockroom staff.					
* For extensive injuries, student is escorted	to the	before	or imme	diately call	
if after					
7. Broken Glassware Box, Dust Pan & Br	room				
* Disposal container for		broken glassware			
* Please use the dust pan and broom to sw	eep up any tiny gl	ass pieces.			
* If it's not,	it doesn't go in he	ere!			
8. Spill Control Center					
* All spills must be					
* () will	neutralize solutior	ns that are acidic /	basic.		
* should be used for a	absorbing spilled s	solvents.			
* For spills larger than a few milliliters, it ma	y be necessary to	evacuate and cal	l x911.		
9. Balance Station					
* Read the BALANCE ETIQUETTE SIGN	* Use piece of _		to weigh solid	ls.	
* to help the t	ransfer.				
* Bring the container you're transferring into	o – DO NOT walk	around the lab witl	h	_	
<u> </u>					
* Always the balance p	oan and			NO SNOW	
STORMS of chemical powder here!					

* After properly disposing of	, wash gla	ssware with	and
, then rinse twice w	ith DI water.		
* This is NOT a waste bin. Only	and	down the dra	ain!
* Note which sinks in this room, if any, ha	ave a flood hose.		
11. Chemical Waste Station			
* Waste bottles are kept here in		(bin to	o catch spills).
* Read the waste label – there may be m	ore than one type o	of liquid waste.	
* Pour <i>into</i> the waste bottle using the	r	provided (not onto the	e bottle. Yes, someone's done
that before. Not cool).			
* if a waste	bottle is full. Don't	let the containers	·
<u>12. Dry Waste Box</u>			
* For solid waste from experiments such	as,	, and _	·
* DO NOT put		in the	dry waste box.
* If it's a liquid, it doesn't go in here!			
* if you are unsure	e of what goes in th	e dry waste box AFT	ER reading the guidelines
above and instructions in lab procedures			
13. Chemical Fume Hoods			
* Minimizes your	to		
* Work with the chemicals at least		into the hood	
* Hood cover/sash should be	to the	or	else!
* DO NOT put in the	hood! * Keep :	surfaces clean – clea	in up spills

14. Reagent Station

* Take only what you nee	ed from bottles	* Keep surfaces clean – clean up spills!	
* Prevent contamination	- DO NOT return .	_to	
* Bring a	and	_ for transfer – DO NOT walk around with a full pipe	t!
* Carefully read labels tw	ice	* Carefully	
<u>15. Disposable Gloves</u>			
* This is a	line of defense	e * Gloves do not make your hands	_!
*	_gloves if you get	chemicals on them	
* Let your TA know if a bo	ox is	also, be a good human and place said empty	box in the
	!		
<u>16. Equipment Room</u>			
* GCs can be hot! Don't le	eave	on top	
* Keep GC/IR kits tidy () when in use	
* Clean up spills when the	ey happen		
* Ask your TA to show yo	u around the instr	rument room.	
17. Rota-vap			
* Used to	samples	s: Round-bottom flask is attached, rotation prevents	boiling over as
is applied	to remove solven	nt, which is collected in the	
* Your TA will		and/or set this up for you the first couple times.	
* Please	when not in us	se. * Be respectful – empty that	_!
18. Water Pump & Wate	<u>r Lines</u>		
* Don't let pumps	, check f	requently	
* Know your in's and out's	s – is w	vater in,is water out	
* Water lines run near		— watch where you point those things!	
* Only use the clamps to	adjust flow – plea	se do NOT	

19. Ring Stands

* Used to secure ______ using ______.

* Stack them in an ______before returning

20. Clamps

* Separate clamps from _____

* Do not leave ______ or other items in this drawer

* If the threads on the clamp become worn and no longer work, please bring it to the

_____ so we can try to fix it

21. Vacuum Tubing

* ______tubing for connecting a vacuum line.

* Ask your TA where to connect to the _____.

* Return tubing when you are finished.

* Do not use for_____; vacuum only.

22. Hot / Stir Plates

* Note that there are separate dials for '_____' and '____' and that different hotplates have

different types of ______ settings.

* Mind where the cord lies to prevent_____

* Set heat on or below ______ setting. These hot plates get ridiculously (unsafely) hot on ____.

* When finished, put them away neatly, no leaning towers of hot plates please.

23 – What is one of your TA's hobbies or interests?

Share something similar with your TA 😊

* 24 & 25 are on the following pages *

26. Lab Coats

* Worn over _____.

* Must be worn with ______ during all experimentation and cleaning.

* Contaminated lab coats are considered waste. Notify your TA and bring the coat to the

_____ if you spill on your lab coat.

* Lab coats are shared with many sections – be considerate and do not leave

____in the pockets!

* Hang up coats neatly on the hanger labeled with the _____.

24. Hazard Definitions

In the case of exposure to any chemical, rinse the affected area immediately for _____and

Copy the following precautions to be taken when handling the following types of chemicals then find the HAZARD TERMS below.



HAZARD TERMS

R S M E B X O A V M C H H W P U W A Y L D Y H J A B T Y K Q SSBFFBMEROJJGUE ΗΟΗDΕΙΑCΧΙSΤRΚΤ SLUGRTIMIOBWOTP NCDCCNYNMYKZSNV FXIOLNFDAYWCAZ F ZBXGYLNRIKLNOTI K O E B C B G X Z R N F P I Т TNVYUWBUOMSRIRB LACHRYMATORTCRW OGKWWATBVANOWIT X D M S E V I S O R R O C U J J R O M A T H L T B J P M I Y ΡΚΟΕΝΗΧСΙJΡΚЈΝΜ

CARCINOGEN	CORROSIVE	FLAMMABLE
HYGROSCOPIC	IRRITANT	LACHRYMATOR
SAFETYFIRST	SLUG	TOXIC

25 – <u>NFPA Labels</u> - Copy the NFPA label description from the bulletin board. Color them in if you have the equipment!

What does NFPA stand for? _____





Classify each of the examples below using the ratings above. Forgive the lack of color, feel free to add your own!



Health Hazard		
Fire Hazard		
Specific Hazard		
Instability		

Part B. Equipment Locker

Clean & Organized Drawer = Safe & Orderly Experiments 😊

- Find the locker **equipment list** in addition to pictures of glassware and the pictures of the drawers.
- Check that all your equipment is present in your drawer.
- There should be no missing, extra, dirty, or broken items after lab.
- At the end of each lab, pick up each item from your drawer to make sure there are no cracks or chemicals present, even if you didn't use that equipment that day (things can fall into your drawer!).

Sketch and name at least three items from your drawer that are new to you

Lab Basics - Balances, Pipets, & Solutions, oh my!

Part C. Making Solutions

- 1. Find the molecular mass, boiling point, melting point, density, and safety hazards in the Material Safety Data Sheet (MSDS) in the lab.
- 2. Enter the chemical name and amount measured for the solid in milligrams and for the liquid in milliliters.
- 3. GO TO THE PROCEDURE. Enter data on the following pages AND below.

	Chemical Name	Amount	Molecular mass (g/mol)	Moles (mmol)	Boiling or melting point (°C)	Density (g/mL)	Safety Hazards
Solid Part C.1		mg					
Liquid Part C.2		mL					

Part C.1. PROCEDURE - Measuring Solids

- The digital balances read in grams (g) so you'll need to convert to milligrams (mg).
- There are 1000 milligrams in every gram $(10^3 \text{ mg} / 1 \text{ g})$ or $(1 \text{ g} / 10^3 \text{ mg})$.
- Convert 50 mg into grams: 50 mg = _____ g

Use the "Balance Etiquette" on next page - weigh approx. 50 mg (± 10 mg) of your assigned solid...

• Note $(\pm 10 \text{ mg})$ means your measurement can be 10 mg higher or lower than 50 mg.

40 mg = _____ g 60 mg = _____ g

Record the exact mass measured: _____ mg

• Convert mass (mg) to millimoles (mmol) using the molecular mass of the solid.

_____ mmol

- Copy relevant data to the table above.
- Transfer the solid into a small beaker or Erlenmeyer flask labeled with chemical name.

PROCEDURE - Balance Etiquette

Line at the balance? Introduce yourself to the person in front of you and they'll tell you when they're done!

- 1. Bring a spatula or scoopula & labeled container into which to transfer the solid.
- 2. Fold the weigh paper in half on the diagonal and place on the balance pan
- 3. Tare (zero) the balance
- 4. Add the solid onto the paper and record the mass
- 5. Place extra solid on separate weigh paper, not back into the original container
- 6. Close the container
- 7. Transfer the solid into the labeled container (ex. Erlenmeyer flask or beaker)
 - DO NOT walk around the lab with solid on weigh paper!
- 8. CLEAN! Brush spilled solids onto separate weigh paper then transfer to the solid waste.
 - Dispose of weigh paper in the trash
 - Wipe down the counter with a wet sponge or paper towel

Part C.2. PROCEDURE - Measuring Liquids and Making Solutions

Use the steps below to measure 15.0 mL of deionized water (diH₂O) and mix with your solid...

- 1. **Label** the 25-mL graduated cylinder from your drawer (diH₂O).
- 2. Locate the diH_2O reagent bottle in the **fume hood**.
- 3. **Pour a little under 15 mL** from diH₂O reagent bottle into the large grad cylinder in your drawer.
- 4. Use a glass pipet and bulb to remove or add smaller volumes of water to get to exactly 15.0 mL.
 - Transfer any excess water into a separate container, NOT back into the reagent bottle.
- 5. Convert volume (mL) to millimoles (mmol) using the density and molecular mass of water.

mmol

- 6. Make the solution: Transfer the diH₂O into the container with the weighed solid from Part C.1.
 - Swirl to dissolve completely.
 - What's happening at a molecular level as the solid seems to disappear?
 - Sketch diagrams of the following forces:

Solid: ionic bonds	Liquid: Hydrogen-bonding	Solution: Ion-dipole interactions

7. Calculate the solution concentration.

Molarity = (mmoles solute) / (mL of solution) = _____ M

Part C.4. PROCEDURE - Mixing Solutions - pluringe and graduated cylinder practice

What's a pluringe and how do I use it? A pluringe is a disposable glass pipet connected to a reusable syringe with rubber tubing (pluringe). This allows the synthetic chemist to obtain small volumes of liquids with reasonable accuracy and precision. The pipet needs only be inserted into the rubber tubing enough to stay put, not all the way, as it becomes difficult to remove later. The pluringe is re-usable as long as it does not come into contact with liquids. It is easy to avoid pluringe contamination:

- DO NOT turn the pluringe upside-down when liquid is in the pipet.
- Obtain only volumes smaller than the pipet or pluringe itself.

Only the glass pipet touches liquid, not the rubber tubing or syringe. Pipets can go in the glass waste if the pipet was only used for water. Broken syringes can be exchanged in your lab – find the little plastic cup with spare good syringes.

Part C.4. PROCEDURE - Mixing Solutions

- Pair up with another student to share solutions.
- Transfer their solution into yours using each of the following tools from your drawer. Record the volume.

Instrument	Volume from instrument (mL)	New volume after mixing (mL)
1-mL pluringe		
- · · ·		
3-mL pluringe		
10-mL grad cylinder		

Group Discussion: What would happen if...

...you quickly pulled back the syringe to obtain liquid? Or if you quickly pushed the syringe to transfer?

...you were transferring a second portion of liquid into a container and dipped the tip of the pipet into the liquid in that container?

...you dispensed more liquid into the pluringe than the capacity of the pipet? What is the volume of a short-stem glass pipet?

...you threw away a glass pipet in the regular trash?

Part C.5. CLEANUP - Keep lab coat, goggles, and mask ON

- 1. Dispose of chemicals in appropriate waste container. Throw away your gloves in the regular trash.
- 2. Wash glassware with soapy water and brush. Rinse with diH₂O. Leave to dry on paper towel.
- 3. Wipe down the counter at your workstation with a wet sponge then dry with a paper towel.
- 4. Ask your TA for a community cleanup task.
- 5. Put away glassware in your drawer, leaving upside-down to dry if necessary and possible!
- 6. Wash your goggles with soapy water. Rinse with diH₂O, blot to dry, then hang on the pegs for goggles.
- 7. Take off your lab coat, empty the pockets, and place on the appropriate hanger by size.
- 8. Wash your hands.
- 9. Double check with your TA that it's ok to leave then say bye to your lab buddies!

<u>Exp's 1 & 2 = Column Chromatography & Acid-Base Extraction for Excedrin Separation</u> Revell, K. D. *Journal of Chemical Education.* **2011**, *88*, 1413.



Functional Groups, Intermolecular Forces (IMFs), & Polarity

- Ion-dipole

- Hydrogen bonding (H-bond)







- Dipole-dipole



<u>CHROMATOGRAPHY</u> →

(S) <u>Stationary Phase</u> = SiO_2

(M) <u>Mobile Phase</u> = Organic Solvents Hexanes, Ethyl Acetate (EtOAc), & Acetone

- Order of Separation

- Degree of Separation

Thin-Layer Chromatography (TLC) Analysis – ACE, ASP, & CAF Standards



TLC Analysis 6h Standardstion

LABORATORY EXPERIMENT



TLC analysis of fractions after column chromatography (Example from Revell article)

Prepare Solvent, then Pack & Load the Column for Liquid Chromatography



Run the Column



9. Combine & concentrate appropriate fractions in vacuo (rotary evaporator, rota-vap)

10. Percent (%) Recoveries from Excedrin tablet (~675 mg)

1 tablet = 250 mg ACE, 250 mg ASP, 65 mg CAF, plus inactive ingredients

Experiment Overview



(1) Extraction / separation

(2) Analysis: TLC & IR

ACIDIC & BASIC FUNCTIONAL GROUPS (FGs)

(accepts or donates protons, H⁺)

- Carboxylic acids, phenols, amines

vs. NEUTRAL FGs

(not acidic or basic)

- Esters, alcohols, hydrocarbons



- Phenols



- Amines



Two-Base Extraction of Excedrin's Active Ingredients



Extraction Flow Chart

- What's where and when?!



LIQUID-LIQUID EXTRACTION

- 2 immiscible liquids have distinct densities and polarities
- Component(s) have a solubility preference for one layer: aqueous (AQ) & organic (ORG)
 - PARTITION COEFFICIENT = ratio of solubilities in 2 solvents

Procedural Lingo...

Basic Extraction

"Extract the ORG layer w/ AQ base" Acidic Extraction "Extract the AQ layer w/ ORG solvent" Washes

"Wash the ORG layer with sat. NaCl (brine)"







1. Isolation of Aspirin



2. Separation of Acetaminophen & Isolation of Caffeine



strong base extraction

3. Isolation of Acetaminophen

- ACE typically does not precipitate upon addition of HCI





Analysis of Excedrin Components via TLC and IR Spectroscopy

- 1. TLC same as Exp 1! TLC mobile phase = 1:2 Hex / EtOAc w/ 1% HAc
 - dilute solids & check for purity before IR analysis

ASP ACE CAF

2. IR analysis = detection of bonds within specific functional groups



- Frequency of vibration depends on bond length

- Longer bonds have greater distance to expand/contract = slower frequency, \Downarrow_v



Antisymmetric stretching



Out-of-plane bending

Identify the functional groups and expected IR ranges from IR table on Canvas.



Conjugated vs. saturated double bonds





CHEM 8M, NMR Spectroscopy for Structure Determination



ethyl propanoate ... says who?!

NMR =

Nuclear Magnetic

Resonance

INSTRUMENTATION

ATOMIC NUCLEI, ex. ¹H atoms



Time (s) FID δ (ppm)

Spectrum

¹H NMR is used determine or confirm a chemical structure

- Chemical Equivalence: How many different types of protons are in the molecule?
- Chemical shift: What is the chemical environment of those protons?
- Integration: How many H's are represented by each signal?
- Splitting: How many H's are nearby?



FIGURE 22.11 ¹H NMR spectrum of ethyl propanoate at 200 MHz.

One type of "NMR-active" nucleus can be observed in an experiment at a time: ¹H, ¹³C, ¹¹B

- ¹H NMR = "proton NMR" is most common / useful because ¹H isotope > 99% abundance.
 - Small amount of sample (< 10mg)
 - \circ Quick run time (< 5 min).

0

- Chemical shift range: 0 12 ppm
- ¹³C NMR = "carbon NMR" is also widely used.
 - ¹³C isotope < 1% abundance
 - more sample (< 30mg)
 - longer run time (30+ min)
 - Chemical shift range: 0 220 ppm



TABLE 21.1 Deuterated solvents for NMR spectroscopy			
Solvent	Structure	Residual ¹ H signal (ppm)	¹³ C chemical shift (ppm)
Chloroform-d	CDCl ₃	7.26 (singlet)	77.0 (triplet)
Acetone-d ₆	$CD_3(C=O)CD_3$	2.04 (quintet)	29.8 (septet), 206.5 (singlet)
Deuterium oxide	D ₂ O	4.6 (broad singlet)	—
Dimethyl sulfoxide-c	$I_6 CD_3(S=O)CD_3$	2.49 (quintet)	39.7 (septet)

CHEMICAL SHIFTS, EQUIVALENCY, & CHEMICAL ENVIRONMENTS

- Nearby atoms / FG's cause the chemical shift to be higher (deshielded) or lower (shielded)
 - Shielding = electron clouds 'protecting' the nucleus from experiencing the magnetic field of the NMR instrument
- Pi bonds, electron-withdrawing groups (EWG's), & electronegative (EN) atoms cause deshielding

Deshielded H's

Strong δ^+ ↑Chemical shift Shielded H's

Weak $\delta^{\scriptscriptstyle +}$

 \downarrow Chemical shift



CH₃OCH₃



Si(CH₃)₄

Relate the terms below to propanoic acid...

- Chemical Equivalence

- Integration

- Chemical Shifts



How to Predict ¹H NMR Spectra



- 1. Chemical Equivalence: How many ¹H NMR absorptions are expected?
 - How many types of non-equivalent protons?
- 2. Integration: What is the ratio of peak areas expected upon integration of the spectrum?
- 3. Approximate the **chemical shift** for each type of proton using the NMR correlation tables.
 - Table 22.2...H-C-("Group") used to assign most appropriate chemical shift range
 - Tables 22.3 and 22.4 are used to calculate an expected chemical shift based on nearby atoms
 - i. Use Table 22.4 for H's connected to a benzene ring; use Table 22.3 for all other H's
- 4. Splitting & (n+1) rule: Number of H's adjacent to each set of H's

H₃CO[~]

0

Proton Splitting aka ¹H-¹H Coupling and the "n+1 Rule"

• Non-equivalent signals (neighboring proton nuclei) effect each other's peak shape.



N = number of hydrogens on adjacent carbons (neighbors)




Frequency (Hz)

¹H spectrum of ethyl propanoate has two 2H Common NMR terminology.

1 two 3H triplets (A&D).



Signal D is a 3H triplet at 1.3 ppm: " δ 1.1 (3H triplet)"

- All D nuclei spin against B field with central chemical shift @ 1.3 ppm
- D has 2 neighbors (signal C, n = 2) that are 50 / 50 with & against B field



FIGURE 22.12 ¹H NMR spectra of *tert*-butyl acetate in the region from 0 to 500 Hz at (a) 60 MHz and (b) 200 MHz. The chemical shift of each signal is the same regardless of the spectrometer frequency.

Preparing for NMR Problem Set

- watch the "How to do the NMR problem set" video (Canvas) and use the template below to follow along with the problem-solving process.

... Predict the number signals, integration, splitting, and chemical shift (range and calculated)

Br	Signal	Integration	Splitting	Chemical Shift Range (ppm)	Calculated Chemical Shift (ppm)*
Br					

Signal	Integration	Splitting	Chemical Shift Range (ppm)	Calculated Chemical Shift (ppm)*

<u>Structural Elucidation</u> - Propose a structure from chemical formula, IR, & ¹H NMR

Watch the "How to do the NMR problem set" video (Canvas) and use the template below to follow the problemsolving process.

- Calculate degrees of unsaturation (total number of pi bonds and/or rings) from the molecular formula
- use IR to identify possible functional group(s),
- draw structure fragments based on signal information (2H doublet = CH₂CH),
- then put the fragments together in the **final structure**.
- Check your work by calculating chemical shifts of the proposed structure.

C₃H₆O; **IR** 1720 cm⁻¹, 2900 cm⁻¹; ¹H **NMR** – δ 2.0 ppm (3H singlet)

C₅H₁₂O; **IR** 3300 cm⁻¹, 2900 cm⁻¹; ¹H **NMR** – δ 4.0 ppm (1H broad singlet), δ 3.5 ppm (2H triplet), δ 1.6 (1H nonet), δ 1.5 ppm (2H quartet), δ 0.9 (6H doublet)

What's the resonance frequency, Kenneth?

See optional/supplemental video (Canvas) on what a chemical shift really is.

Chemical shift (ppm) = <u>frequency of proton relative to TMS</u> Maximum operating frequency of instrument (MHz)

Ex. *tert*-butyl acetate is dissolved in a solution of CDCl₃ (w/ 0.1 % TMS) and ¹H NMR spectra are run on two different instruments.

(a)



FIGURE 22.12 ¹H NMR spectra of *tert*-butyl acetate in the region from 0 to 500 Hz at (a) 60 MHz and (b) 200 MHz. The chemical shift of each signal is the same regardless of the spectrometer frequency.

Binder

¹H NMR Introductory Problems

1. How many ¹H NMR absorptions are expected for each compound? In other words, how many non-equivalent protons are in each compound? *Pro-tip: Draw the skeletal structure then add H's.*

(a) CH_3CH_2CI (b) $(CH_3)_2CHOCH_3$ (c) $NO_2CH_2CH_2CH_3$

(d) Benzene (e) 2-Methyl-1-butene (f) *trans*-3-hexene

2. Each compound below has only 1 peak in its ¹H NMR spectrum. Use the NMR tables to predict the chemical shift range (ppm) where each compound should absorb.



3. Consider the structure of *p*-xylene below to answer the following.



- (a) How many peaks (absorptions) should *p*-xylene have in its ¹H NMR spectrum?
- (b) What ratio of peak areas would you expect on integration of the spectrum?
- (c) Use the NMR tables to approximate the chemical shifts of the signals in *p*-xylene.

4. Summarize in your own words *how* to find the number of signals, integration, and chemical shift(s) in any given molecule.

5. Draw (invent) or look up a molecule you know by name online (ex. medicine). Draw that structure below. Predict the number of ¹H NMR signals, integration, and chemical shifts of that molecule! Share your findings to your labmates.

NMR Problem Set					
Pages 1-4	Pages 5-8	Page 9			
#1 Predict spectra from given	#2-4 Structure elucidation from	#5 Summary / Overview			
structure	given spectra				

- **1.** Predict the ¹H NMR spectrum for compounds (b) (f) below.
 - Add in the missing hydrogens
 - Determine the **number of** ¹H NMR signals in each compound below.
 - Label each **type of proton** (A, B, etc.) and fill in the tables provided on the following page with a row for each type of proton (signal).
 - Predict the **integration** (#H's), the **chemical shift range** (ppm), and **calculated chemical shift** (ppm), that you would expect for each signal.
 - Determine the **splitting pattern** for each signal using the <u>*n*+1 rule</u>, where *n* is the number of H's on adjacent carbon atoms.
 - Splitting patterns: singlet, doublet, triplet, quartet, pentet, sextet, heptet, or multiplet.

= 2.3 ppm

- There is no splitting through heteroatoms (O, N, S, etc.) all OH's are singlets!
- Sketch the predicted ¹H NMR spectrum of each compound and use <u>nmrdb.org</u> to check your work
 - Approximate the central chemical shift of the signal on the x-axis and *draw* each splitting pattern.
 - Don't worry about incorporating integration (peak size).



(1a) Worked example: Watch the "How to do the NMR Problem Set" video on Canvas

Br_CH ₃	Signal	Integration	Splitting n+1 rule	Chemical Shift Range (ppm)	Calculated Chemical Shift (ppm)*
A ^H Br	A	1H	quartet	2.1 – 4.5	5.7
	В	ЗH	doublet	0.8 - 1.9	2.3

Calculate expected chemical shifts using *Table 22.3 of NMR tables – alkyl H's

Calculation of Signal A		
Methine (CH) +	2 x Br (alpha substituent)	
1.5 ppm +	(2 x 2.1)	= 5.7 ppm
Calculation of Signal B		
Methyl (CH ₃) +	2 x Br (beta substituent)	

+ (2×0.7)

0.9 ppm

p1

0

Tables for #1 – create as many rows as the number of signals in each compound, read page 1 for more!

(46)	Signal	Integration	Splitting	Chemical Shift Range (ppm)	Calculated Chemical Shift (ppm)*

Sketch of ¹H NMR spectrum:

12

Chemical Shift (ppm)

Sketch of ¹H NMR spectrum:

Tables for #1 – create as many rows as the number of signals in each compound, read page 1 for more!

(1d)	Signal	Integration	Splitting	Chemical Shift Range (ppm)	Calculated Chemical Shift (ppm)*

Note: Use correlation Table 22.4 for aromatic H's and Table 22.3 for alkyl H's

Sketch of ¹H NMR spectrum:

12

Chemical Shift (ppm)

0

(1e)	Signal	Integration	Splitting	Chemical Shift Range (ppm)	Calculated Chemical Shift (ppm)*
, i dy					

Note: Use correlation Table 22.4 for aromatic H's; no calculation for OH

Sketch of ¹H NMR spectrum:

Tables for #1 – create as many rows as the number of signals in each compound, read page 1 for more!

	Signal	Integration	Splitting	Chemical Shift Range (ppm)	Calculated Chemical Shift (ppm)*
(1f)					
о о о он					

Sketch of ¹H NMR spectrum:

- 2. Structural Elucidation: Draw structures for compounds that meet the following descriptions.
 - 1. Calculate **degrees of unsaturation** from the **molecular formula** to determine the number of pi bonds or rings
 - 2. Use IR to identify potential bonds & functional group(s)
 - 3. Draw structure fragments for each ¹H NMR signal using the integration & splitting
 - a. Ex. **2H doublet** = C**H**₂CH, where the bold H's are determined by integration (2H) and the neighboring CH group is determined by splitting (doublet)
 - b. Incorporate functional groups into fragments where possible using chemical shifts
 - i. Ex. δ 4.0 (2H doublet) = OCH₂CH
 - ii. the 2H's are deshielded by the electronegative oxygen, resulting in a higher chemical shift
 - 4. Put the **fragments** together in the final structure.

Formula	Degrees Unsat.	Suspected bonds / FG(s) from IR	Fragment(s) from ¹ H NMR	Final Structure
C3H6O	1	1720 cm ⁻¹ C=O 2900 cm ⁻¹ Alkane sp ³ C-H	δ 2.0 ppm (6H singlet) 6H's are next to 0 H's CH ₃ C _(no H's) Symmetry - 2 x CH ₃ groups	o
C3H7CI	0	2900 cm ⁻¹ Alkane 800 cm ⁻¹ C-Cl bond	δ 3.7 (6H doublet) ^ next to Cl Doublet = 6H's next to one H (CH ₃) ₂ CHCl δ 1.4 (1H septet) ^ farther from Cl Septet = 1H next to 6 H's (CH ₃) ₂ CHCl	(C <mark>H₃</mark>)₂C <mark>H</mark> CI

Worked examples: Watch the "How to do the NMR Problem Set" video on Canvas

Formula	Degrees Unsat.	Suspected bonds / FG(s) from IR	Fragment(s) from ¹ H NMR	Final Structure
			δ 4.1 (2H quartet)	
(2a)		1760 cm ⁻¹	δ 2.3 (2H quartet)	
$C_5H_{10}O_2$		2900 cm ⁻¹	δ 1.3 (3H triplet)	
			δ 1.1 (3H triplet)	
			δ 4.0 (2H quartet)	
(2b)		1760 cm ⁻¹	δ 2.0 (3H singlet)	
C4H8O2		2900 cm ⁻¹	δ 1.4 (3H triplet)	
			δ 3.8 (1H septet)	
(2c)		3300 cm ⁻¹ 2900 cm ⁻¹	δ 2.0 (1H broad singlet)	
031180			δ 1.2 (6H doublet)	

Formula	Degrees Unsat.	Suspected bonds / FG(s) from IR	Fragment(s) from ¹ H NMR	Final Structure
			δ 7.85 (2H doublet)	
		3050 cm ⁻¹	δ 7.28 (2H doublet)	
(2d) Bonus Problem:		1675 cm⁻¹	δ 2.9 (2H quartet)	
C ₁₀ H ₁₂ O		1620 cm ⁻¹	δ 2.0 (3H singlet) *benzylic	
		850 cm ⁻¹	δ 1.5 (3H triplet)	

3. Determine the **structure** of a compound with molecular formula $C_4H_{10}O$ and the integrated ¹H NMR spectrum below (**rel. area** gives the *ratio* of H's per signal). Identify the splitting pattern from the expanded shapes in the spectrum (**m = n + 1**, where **m** = multiplicity, triplet or quartet and **n** = number of neighboring H's).

Show your work: draw fragments of the molecule in addition to the full structure.
 ↓ This box gives the relative area (integration or ratio of H's) of both peaks, identified by chemical shift.



4. Propose a structure for a compound with formula $C_5H_{10}O_2$ and a sharp IR absorbance at 1735 cm⁻¹. Identify the splitting pattern from the expanded shapes in the spectrum (**m** = **n** + **1**, where **m** = multiplicity, triplet or quartet and **n** = number of neighboring H's). The signal at 2.0 is a singlet.

Chem. Rel. shift area 1.22 6.00 2.01 3.00 4.99 1.00 Intensity 10 8 7 3 2 9 6 5 4 0 ppm 1 Chemical shift (δ)

•

- Show your work: draw fragments of the molecule in addition to the full structure.
- \downarrow This box gives the relative area (integration) of all peaks, identified by chemical shift.

(a) Predicting Spectra

(b) Structural Elucidation

CHEM 8M, Experiment 3 – Oxidation of Benzhydrol

- Reactions using Phase Transfer Catalysts (PTC)
- TLC, IR, and ¹H NMR Analysis of Benzhydrol & Benzophenone

Oxidation of Benzhydrol with Bleach using Phase Transfer Catalyst



Phase Transfer Catalysis (PTC)





Monitoring Reaction Progress by TLC



Reaction Work-up: remove AQ, wash with brine then water, dry, filter, rota-vap

IR Analysis

 $\begin{array}{c} \mathsf{OH} & \longrightarrow & \mathsf{O} \\ \mathsf{Ph} & \longrightarrow & \mathsf{Ph} & \overset{\mathsf{O}}{\vdash} \mathsf{Ph} \end{array}$

¹H NMR Analysis = assign each set of protons on structure to signal on spectrum

Predict spectrum

1. Look for symmetry – equivalent protons and for asymmetry – non-equivalent protons





2. Integration: How many of each type of proton?

3. Identify chemical shift ranges in benzhydrol & benzophenone



Benzhydrol



Benzophenone

4. Calculate expected chemical shifts using chemical shift correlation tables or online predictor tool

 Table 3. ¹H NMR Analysis of Benzhydrol

с _{ННО Н} Е	Signal	Integration (# of H's)	Expected Chemical Shift (ppm)	Observed Chemical Shift (ppm)
"H	Α			
A _H	В			
Development	С			
Benzhydrol	D			
	E			

Correlate / assign to signals on given spectrum - Integration lines = Curves on or above peaks, height = relative ratio of H's



Table 4. ¹H NMR Analysis of Benzophenone

C _H O B _H ↓ ↓ ∧	Signal	Integration (# of H's)	Expected Chemical Shift (ppm)	Observed Chemical Shift (ppm)
	Α'			
AH, A	B'			
Benzophenone	C'			

Correlate / assign to signals on given spectrum



Supplemental Reading

Klein, D. "Organic Chemistry, 2nd Edition"; **Mohrig**, J. R.; *et. al.* "Techniques in Organic Chemistry, 4th Edition." See also Exp 3 pre-lab videos on Canvas

Oxidation reactions	Klein 12.10
¹ H NMR	Klein 15.1-6 or Mohrig Chapter 22.1-22.7
Extraction	Mohrig Chapter 10
TLC	Mohrig Chapter 18

F. Experimental Methods Writing Worksheet - provided in lab ©

1. Draw the **reaction scheme** by hand (no copy/paste) and list the **name of the product**. *The reaction scheme includes reactant, reagents over arrow, solvent under arrow, and product.*

2. What glassware and equipment was used for this reaction (aside from chemicals)?

3. How much **benzhydrol** was used? Convert **mass** to **mmol** (**xx g**, **xx mmol**). Show your work, including units with every value. Calculate or look up the molecular weight of benzhydrol (g/mol) = (mg/mmol).

4. How much bleach (**NaClO**) was used and what was the **concentration** (_____**M**, ____**mL**)? Fill in the blanks and calculate the quantity of **bleach in mmoles.** Show your work. *Recall Molarity* = (*moles / Liter*) ... M = (mol / L) = (mmol / mL).

5. How much *tert*-butylammonium hydrogen sulfate (**Bu**₄**NHSO**₄) was used (**xx g**)? This is a catalyst – include only mass not mmol.

6. Determine the **limiting reagent** then calculate the **theoretical yield** (mmol and mg). Show your work, including units with every value. *Determine the mole ratio in the reaction (x mol benzyhydrol / x mol benzophenone).* Calculate or look up the molecular weight of benzhydrophenone (g/mol) = (mg/mmol).

F. Experimental Methods Writing Worksheet (cont'd)

7. What **solvent** was used in the oxidation reaction and in what **volume**?

8. What was the reaction temperature and time? Was the reaction stirred, refluxed, or standing?

9. What technique was used to monitor reaction progress? What solvent(s) were used during this analysis?

10. List the **identity** and **quantities** of the **chemicals (xx mL)** used in the **reaction work-up**. *Note: quantity of drying agent need not be included.*

11. What additional processes were involved in the final isolation of product?

12. What is the yield of **benzophenone** (_____ **g**, _____ **mmol**, _____ **% yield**)? Fill in the blanks and show your work below, including units on every value.

(a) Convert benzophenone mass (300 mg) to mmol using molecular weight (g/mol) = (mg/mmol).

(b) Calculate percent (%) yield using 300 mg as the actual yield and the th. yield from #6.

% yield = <u>actual yield (mg)</u> x 100% Theoretical yield (mg)

CHEM 8M, Experiment 4 – Synthesis of Fruity Fragrances

Fischer Esterification

- = Reaction of <u>carboxylic acid</u> and <u>alcohol</u> under <u>acidic, refluxing</u> conditions
- "Neat" reaction!





Le Chatelier's Principle – brainstorm!



Hydroxamic Acid Test for Esters



¹H NMR Banana Oil & Apple Oil

- How many different types of protons (signals) do you see?
- What's the splitting pattern for each signal?

 $\sim \sim \sim$

Predicting Relative Chemical Shifts

- What is the expected chemical shift range & calculated value for each set of H's?
- Deshielding increases chemical shifts
 - Electron withdrawing groups (EWGs) or electronegative (EN) atoms
 - o Proximity effects





Assigning Banana Oil's Spectrum Signals A-E represent sets of protons in the product



Signal	Integration (# of H's)	Expected Chemical Shift (ppm)	Observed Chemical Shift (ppm)	# of 3-bond Neighbors (n)	Multiplicity/ Splitting (m = n+1) ex. Singlet, doublet, etc.
A			4.0		
В			2.0		
С			1.6		
D			1.4		
E			0.9		

¹H NMR Analysis of Hexyl Acetate (Apple Oil)

Signal	Integration (# of H's)	Expected Chemical Shift (ppm)	Observed Chemical Shift (ppm)	# of 3-bond Neighbors (n)	Multiplicity/ Splitting (m = n+1)
Α	2H		4.0		triplet
В	3H		2.0		singlet
С	2H		1.5		pentet
D	6H		1.2	(Many)	multiplet
E	3H		0.8		triplet





CHEM 8M, Experiment 5 – Synthesis of Aspirin

- Reaction Set Up, Work Up, Chemical Tests, IR,
- ¹H NMR Analysis and intro to ¹³C NMR



Reaction Work Up

1. Cool for ~1 min

2. Add water to quench in warm water bath, 5 min



3. Crystallize

- (a) Cool to RT, transfer to beaker
- (b) Cool in ice bath, then scratch inside bottom of beaker, then wait!
- (c) No crystals after ~5 min? Add seed crystal
 - Wait at least 5 min after adding seed crystal to filter

4. Vacuum Filtration



Ferric Chloride Test for Phenols



IR Spectra of Salicylic Acid & Aspirin



¹H NMR Analysis

Resonance effects in aromatic rings: Predicting relative chemical shifts without calculations

- Deshielding increases chemical shifts



- Shielding decreases chemical shifts

EDG shields the ortho & para H's





*11.5 ppm broad singlet expected by not observed in this particular spectrum

What differences are expected / observed in the ¹H NMR of aspirin & its precursor?





¹³C NMR

- Exploring carbon nuclei of ¹³C isotopes (< 1.1% abundance)
 longer experiment, requires more sample
- ¹³C nuclei resonate at higher chemical shifts (10-220ppm) than ¹H nuclei (0-12ppm)
- Similar deshielding effects to ¹H NMR



Acetylsalicylic Acid (Aspirin)





0	Chemical Shift (Observed ppm)	Assignment(s) (A' – J')	Expected Chemical Shift Range (ppm)
	169 & 170		
C J OH	152		
D'F' E'	125 – 135 (4 peaks)		
Acetylsalicylic Acid (Aspirin)	122		
	20		

CHEM 8M, Experiment 6 – Colorful Chemistry

Synthesis & Application of	Azo Dyes	
 Nature of color 	- Dye to Fabric Interactio	ns - Diazonium Coupling
Outcomes: Observe effect	s of dye structure, fibers, and	d metals (mordants) on appearance

UV-Visible Spectrum

Purple	Blue	Green	Yellow	Orange		Red
400	450	500	550	500	650	700
λ(nm)						

What is it about the structure of the dye that causes it to appear (emit) a specific color?



Fabric Fibers - Polymers with repeating units of...



Mordant Dyeing

- Fabric strip is pre-treated with coordinating metal: $\rm Cu^{2+},\, Al^{3+},\, or\, Fe^{2+}$
- Pre-made swatches with copper (II) sulfate, aluminum potassium sulfate, or iron (II) sulfate
- Use same ingrain dye procedure above with 'mordant fabrics'
- How does this effect dye fiber interactions?



Direct Dyeing



Part D.1. Synthesis



Part D.2. Vat Dyeing with Indigo


Vibration	Position (cm ⁻¹)	Intensity*	Notes
Alkanes C-H stretch	2990 – 2850	m to s	
Alkenes =C-H stretch C=C stretch =C-H bend	3100 – 3000 1680 – 1620 (sat.) 1650 – 1600 (conj.) 995 – 685	m w to m s	See Table 2 for detail
		-	
alkynes ≡C-H stretch C≡C stretch	3310 – 3200 2250 – 2100	s m to w	
Aromatic Compounds			
C-H stretch C=C stretch C-H bend	3100 – 3000 1625 – 1440 900 – 680	m to w m to w s	Hidden in fingerprint region See Table 2 for detail
Alcohols** O-H stretch	3550 – 3200	br, s	Hydrogen bonded (typical)
Amines N-H stretch	3550 – 3250	br, m	Primary (two bands) Secondary (one band)
Nitriles C≡N stretch	2280 – 2200	S	
Aldehydes C-H stretch	2900 – 2800 & 2800 – 2700	S	H-C=O Fermi doublet
C=O stretch	1740 – 1720 (sat.) 1715 – 1680 (coni.)	S	
Ketones C=O stretch	1750 – 1705 (sat.) 1700 – 1665 (conj.)	s	
Esters** C=O stretch	1765 – 1735 (sat.) 1730 – 1715 (conj.)	S	
Carboxylic Acids** O-H stretch C=O stretch	3200 – 2500 1725 – 1700 (sat.) 1715 – 1680 (conj.)	br, m to w s	
Amides N-H stretch	3500 – 3150	m	Primary (two bands)
C=O stretch	1700 – 1630	S	Secondary (one band)

 Table 1. Characteristic IR Absorption Peaks of Functional Groups^{*}

Table 1 cont'd			
Vibration	Position (cm ⁻¹)	Intensity	Notes
Anhydrides**			
C=O stretch	1850 – 1800 & 1790 – 1740	S	
Acid Chlorides			
C=O stretch	1815 – 1770	S	
Nitro Compounds			
NO ₂ stretch	1570 – 1490 & 1390 – 1300	S	
Thiols [†]			
R-S-H stretch	2550 – 2600		
Alkyl & Aryl Halides [†]			
C-F stretch	1000 – 1400		Hidden in fingerprint region
C-CI stretch	< 600 – 840		
C-Br stretch	< 700		
C-I stretch	< 600		

* Abbreviations: s = strong; m = medium; w = weak; br = broad; sat. = saturated; conj. = conjugated ** Alcohols, Esters, Carboxylic Acids, and Anhydrides also absorb in the fingerprint region due to the C-O stretch (1300 – 1000, s).

Table 2. Out-of-Plane C-H Bending Vibrations in Alkenes and Aromatics

Alkene Structure	Position (cm ⁻¹)	Phenyl Structure	Position (cm ⁻¹)
Mono-substituted R H >=	997 – 985 & 915 – 905	Mono-substituted	770 – 730 & 720 – 680
Disubstituted, <i>trans</i> R H >=-< H R	980 – 960	Disubstituted, <i>ortho</i>	770 – 735
Disubstituted, <i>cis</i> R − R ≻=– H − H	730 – 665	Disubstituted, <i>meta</i>	810 – 750 & 725 – 680
Disubstituted, symm. $R \xrightarrow{H}$ $R \xrightarrow{H}$ $R \xrightarrow{H}$	895 – 885	Disubstituted, <i>para</i>	860 800
Trisubstituted R R R H	840 – 790	R-	800 – 800

^{*} Adapted from...Mohrig, J. R.; Hammond, C. N.; Schatz, P. F. "Infrared Spectroscopy" in *Techniques in Organic Chemistry*. Freeman: New York, 2006.

[†] Palleros, D. R. "Infrared Spectroscopy" in *Experimental Organic Chemistry*. Wiley: New York, 2000. p. 688.



FIGURE 22.13 Approximate regions of chemical shifts for different types of protons in organic compounds.

TABLE 22.2	Characteristic ¹ H NMR chemical shifts in C	
-------------------	--	--

Compound	Chemical shift (δ, ppm)
TMS	0.0
Alkanes (C–C–H)	0.8-1.9
Amines (C–N–H)	0.6-4.5
Alcohols (C–O–H)	0.5-5.0
Alkenes ^a (C=C-C-H)	1.5-2.6
Alkynes (C≡C− H)	1.7-3.1
Carbonyl compounds (O=C-C-H)	1.9–3.3
Halides (X–C–H)	2.1-4.5
Aromatic compounds ^b (Ar–C–H)	2.2-3.0
Alcohols, esters, ethers (O–C–H)	3.2-5.3
Alkenes (C=C-H)	4.5-8.5
Phenols (Ar–O–H)	4.0-8.0
Amides (O=C-N-H)	5.5-9.5
Aromatic compounds (Ar–H)	6.5–9.0
Aldehydes (O=C-H)	9.5-10.5
Carboxylic acids (O=C-O-H)	9.7–12.5

a. Allylic protons.

b. Benzylic protons.

Figures from Mohrig's *Techniques in Organic Chemistry*, 3rd Edition.

IADLE 22.4	chemical shifts of aromatic protons in CDCl ₃				
	Base value	7.36 ppm ^a			
Group	ortho	meta	para		
-CH ₃	-0.18	-0.11	-0.21		
$-CH(CH_3)_2$	-0.14	-0.08	-0.20		
-CH ₂ Cl	0.02	-0.01	-0.04		
-CH=CH,	0.04	-0.04	-0.12		
—CH=CHĀr	0.14	-0.02	-0.11		
-CH=CHCO ₂ H	0.19	0.04	0.05		
—CH=CH(C=O)Ar	0.28	0.06	0.05		
Group	ortho	meta	para		
—Ar	0.23	0.07	-0.02		
—(C=O)H	0.53	0.18	0.28		
(C=O)R	0.60	0.10	0.20		
(C=O)Ar	0.45	0.12	0.23		
(C=O)CH=CHAr	0.67	0.14	0.21		
$-(C=O)OCH_3$	0.68	0.08	0.19		
$-(C=O)OCH_2CH_3$	0.69	0.06	0.17		
—(C=O)OH	0.77	0.11	0.25		
(C=O)Cl	0.76	0.16	0.33		
$-(C=O)NH_2$	0.46	0.09	0.17		
—C≡N	0.29	0.12	0.25		
—F	-0.32	-0.05	-0.25		
—Cl	-0.02	-0.07	-0.13		
—Br	0.13	-0.13	-0.08		
—OH	-0.53	-0.14	-0.43		
—OR	-0.45	-0.07	-0.41		
—OAr	-0.36	-0.04	-0.28		
-O(C=O)R	-0.27	0.02	-0.13		
O(C=O)Ar	-0.14	0.07	-0.09		
$-NH_2$	-0.71	-0.22	-0.62		
$-N(CH_3)_2$	-0.68	-0.15	-0.73		
-NH(C=O)R	0.14	-0.07	-0.27		
$-NO_2$	0.87	0.20	0.35		

TABLE 2.2.4 Additive parameters for predicting NMR

a. Base value is the measured chemical shift of benzene in CDCl_3 (1% solution).

Base values					
	Methyl Methylene Methine	0.9 ppm 1.2 ppm 1.5 ppm			
Group (Y)	Alpha (α) substituent	Beta (β) substituent	Gamma (y) substituent		
	$\mathbf{H} - \mathbf{C} - \mathbf{Y}$	$\mathbf{H} - \mathbf{C} - \mathbf{C} - \mathbf{Y}$	$\mathbf{H} - \begin{array}{c} & & \\ \mathbf{C} - \begin{array}{c} \mathbf{C} \\ \mathbf{C} \\ \mathbf{C} \end{array} \\ \mathbf{C} \\ \mathbf{C}$		
—R	0.0	0.0	0.0		
—C=C	0.8	0.2	0.1		
C=CAr ^b	0.9	0.1	0.0		
-C=C(C=O)OR	1.0	0.3	0.1		
—C≡C−R	0.9	0.3	0.1		
—C≡C−Ar	1.2	0.4	0.2		
—Ar	1.4	0.4	0.1		
(C=O)OH	1.1	0.3	0.1		
(C=O)OR	1.1	0.3	0.1		
(C=O)H	1.1	0.4	0.1		
(C=O)R	1.2	0.3	0.0		
(C=O)Ar	1.7	0.3	0.1		
$-(C=O)NH_2$	1.0	0.3	0.1		
(C=O)Cl	1.8	0.4	0.1		
—C≡N	1.1	0.4	0.2		
—Br	2.1	0.7	0.2		
—Cl	2.2	0.5	0.2		
—OH	2.3	0.3	0.1		
—OR	2.1	0.3	0.1		
—OAr	2.8	0.5	0.3		
O(C=O)R	2.8	0.5	0.1		
O(C=O)Ar	3.1	0.5	0.2		
$-NH_2$	1.5	0.2	0.1		
—NH(C=O)R	2.1	0.3	0.1		
NH(C=O)Ar	2.3	0.4	0.1		

TABLE 22.3 Additive parameters for predicting NMR chemical shifts of alkyl protons in $CDCl_3^a$

a. There may be differences of 0.1-0.5 ppm in the chemical shift values calculated from this table and those measured from individual spectra.

b. Ar = aromatic group.

CHEM 108M, Nuclear Magnetic Resonance (NMR)



TABLE 23.1 Characteristic ¹³C NMR chemical shifts in CDCl₃

Compound	Chemical shift (ppm)
TMS	0.0
$\mathbf{CDCl}_{3}(t)$	77
Alkane (C–CH ₃)	7–30
Alkane $(C-CH_2)$	15–40
Alkane (C–CH) and (C–C)	15–40
Carboxylic acids, esters, and amides (C-C=O)	20–35
Allyl (C –C=C)	20–35
Arene (C -Ar)	20–45
Ketones, aldehydes (C -C=O)	30–45
Amines (C–N)	30–65
lodides (C-1)	-5-45
Bromides (C–Br)	25–65
Chlorides (C–Cl)	35–70
Fluorides (C–F)	80–95
Alcohols (C–OH), ethers (C–OR), esters (C–O[C=O]R)	55-80
Alkyne (C= C)	70-100
Alkene (C= C)	110–150
Aromatic	110–160
Nitriles (C≡ N)	110–125
Carboxylic acids, esters, and amides (C= O)	160–180
Aldehydes (C=O)	185–210
Ketones (C–O)	190–220

Figures from Mohrig's *Techniques in Organic Chemistry*, 3rd Edition.

FUNCTIONAL GROUP*	CLASSIFICATION	EXAMPLE	CHAPTER	FUNCTIONAL GROUP*	CLASSIFICATION	EXAMPLE	CHAPTER
R−X∷ (X=Cl, Br, or I)	Alkyl halide	Öl: n-Propyl chloride	7	ю. В В В	Ketone	·O· L 2-Butanone	19
R R R R	Alkene	1-Butene	7, 8	ю. R Н	Aldehyde	·O· H Butanal	19
R—C≡C—R	Alkyne	1-Butyne	9	;;; В , , , , , , , , , , , , , , , , , , ,	Carboxylic acid	·o· H Pentanoic acid	20
R—ÖH	Alcohol	Butanol	12	r ₩ ਸ਼ ₩ X:	Acyl halide		20
R−Ö−R	Ether	Diethyl ether	13	R, O, R	Anhydride	ioi ioi .o. Acetic anhydride	20
R− <u>ÿ</u> H	Thiol		13	°Ö` R ^{⊥⊥} .Q. [−] R	Ester	·o· .o. Ethyl acetate	20
R−S⊤R	Sulfide	<u>S.</u> Diethyl sulfide	13	R R R R R	Amide	UH _z Butanamide	20
\bigcirc	Aromatic (or arene)	Methylbenzene	17, 18	R I R ^{/N} R	Amine	H N Diethylamine	22
* The "R" refers to the	remainder of the co	ompound, usually carbor	n and hydrogen	atoms.			

TABLE 2.1 EXAMPLES OF COMMON FUNCTIONAL GROUPS

Klein, D. (2019) Organic Chemistry, 3rd edition.