Experiment 4 – Preparation of Fruity Fragrances

Learning Objectives

- Perform and understand Fischer esterification reactions
- Apply acid-base extraction in the reaction work-up
- Critical analysis of liquid-liquid extraction technique
- Observe and interpret hydroxamic acid tests for esters
- Interpret infrared (IR) spectra of starting materials and product to determine reaction success
- Predict and interpret ¹H NMR spectra of synthetic banana and apple oils

* Please find "How to Prepare for Lab & Assignments" after the procedure in this doc.

Background: Esterification Reactions

Esters encompass a large family of organic compounds with broad applications in medicine, biology, and industry. Esters are represented by the structure R(C=O)OR', in which R and R' are alkyl or aryl groups. Esters are widespread in nature, occurring naturally in plants and animals. Small esters, in combination with other volatile compounds, produce the pleasant aroma of fruits. A symphony of chemicals is typically responsible for specific fruity fragrances, however, often one single compound plays the leading role. For example, artificial pineapple flavor contains more than twenty ingredients but ethyl butyrate is the major component. Examples of ester flavors and fragrances are shown in **Figure 1**. In contrast to previous experiments where students isolated compounds from plants, in this experiment, students will synthesize these compounds in the lab.



Figure 1. Examples of esters found in essential fruit oils.

Esters are carboxylic acid derivatives commonly synthesized by Fischer esterification (**Figure 2**). A carboxylic acid is reacted with an alcohol in the presence of catalytic amounts of mineral acids such as sulfuric or hydrochloric acids under refluxing conditions (heat to boiling). This reaction is reversible and thus is limited by its equilibrium constant and dictated by Le Chatelier's Principle. A large excess of either reactant pushes the equilibrium to favor products, thus increasing the yield. Constantly removing the product will also increase the yield, though this is not always possible or practical. A Fischer esterification is only recommended with primary and secondary alcohols and unhindered carboxylic acids. Steric hindrance near the reaction center slows down the esterification.

Figure 2. General scheme for a Fischer esterification reaction.

In this experiment, students prepare either banana oil from acetic acid and isoamyl alcohol or prepare sour apple oil from acetic acid and *n*-hexanol (**Figure 3**). Incidentally, isoamyl acetate is also the alarm pheromone of the honeybee and thus, it should be kept away from beehives! The reaction is performed at the microscale level using a Fischer esterification under refluxing conditions. A round-bottom flask is topped with a water-cooled condenser. The contents of the flask are heated to the boiling. Vapors travel up and inside the reflux condenser, where they condense back to a liquid and fall back into the round-bottom flask. This allows the system to remain open while heating, but without losing any reaction components. This particular reaction is run "neat" or without solvent. The acid and alcohol are both liquids and act as the solvent. It is important that the reaction flask does not run dry and that *cold water* be running through the condenser at all times.



Figure 3. Reaction schemes for fruity fragrance synthesis via Fischer esterification

An acid-base extraction is performed with an aqueous bicarbonate (baking soda, NaHCO₃) solution to separate the ester from the unreacted acetic and sulfuric acid. This weakly basic solution also contains NaCl to improves phase separation in liquid-liquid extractions. This NaHCO₃-NaCl solution has a high ionic strength and draws residual water out of the organic layer. No additional organic solvent is necessary for the acid-base extraction because these esters are liquids and separate from the aqueous solution as an immiscible layer.

Time and quantity permitting, the ester product may be purified by microscale distillation using a Hickman still and a water-cooled condenser. Acid-base extraction is not an applicable method for separation of alcohol from ester, neither of which are acidic or basic! Column chromatography would be effective for separation since the alcohol and ester have very different polarities. The product is analyzed by IR and the hydroxamic acid test for esters and with comparison to alcohol starting materials. ¹H NMR spectra of banana and apple oil are provided for analysis.

PROCEDURE

Procedure Diagrams must be complete in your lab notebook before you can start the lab (see worksheet) The statements in quotes are provided to give you guidance in writing the experimental methods section. One well-written sentence can explain an entire paragraph's worth of information!

1. Reaction Preparation and Set-Up

"To a 15-mL RBF was added...[chemical names (mmol, mL)]...and heated to reflux for 1 hour."

Pre-heat a sand bath on a hot plate at a medium setting. You may set this up as soon as you enter the lab, before the TA's pre-lab talk. Dispense 10 mmoles of the desired alcohol (isoamyl alcohol or *n*-hexanol) and 40 mmoles of glacial acetic acid into a 15-mL round-bottom flask (RBF) using a glass pipet and pluringe. Convert the *mmole quantities into volume (mL) before lab.* Add 3 drops of sulfuric acid and magnetic stir bar then attach a microscale water-jacketed condenser (figure on next page). Be sure the water is running through the condenser and reaction is stirring before heating. Heat to reflux with stirring in the sand bath and allow the reaction to reflux for one hour.

2. Reaction Work-up

"The reaction was quenched and washed with..."

Carefully lift the apparatus from the heat and allow the mixture to cool to ambient temperature. Disassemble the apparatus and turn off the water (clamp), but keep the hoses attached. Do not wash the condenser at this stage, as it may be used later. Transfer the liquid to a 16 x 125 screw-cap test tube with a pipet. Rinse the RBF with 2 mL of 5% NaHCO₃ in 15% NaCl solution. Slowly transfer the rinse to the screw-cap test tube. Stir the mixture with a microspatula until gas evolution (carbon dioxide) has subsided. Cap the tube and invert it several times to mix the layers. Frequently vent the system to release the pressure by momentarily unscrewing the cap. Let the system settle for about 10 minutes.

Use a pipet to transfer the lower aqueous layer to a labeled test tube. Keep this until the end of the experiment then discard it. Wash the organic layer remaining in the test tube twice with 1 mL of the NaHCO₃-NaCl solution. Invert and vent well in each wash. Collect the aqueous washes in the same labeled test tube as before.

"The combined organic layers were dried (Na₂SO₄) and filtered to afford..."

Remove any visible water from the product with a pipet. Dry the organic layer by adding a small microspatula-ful of anhydrous Na₂SO₄. Note that this drying agent is more granular than MgSO₄ and will create a similar but not identical snow-globe effect when sufficient drying agent is added. It may be necessary to add more, however, this may affect how much liquid can be obtained after filtration. Allow the product to dry for 5 minutes with occasional swirling. Filter using a pipet loosely packed with a small piece of cotton into a pre-weighed, labeled vial to obtain the mass of the crude product.

(No purification by distillation)

CHEM 8M, Binder 4. Analysis: Hydroxamic Acid Test

"Product formation was confirmed (or not) by the hydroxamic acid test for esters."

Perform this test with starting materials (alcohol and acetic acid), ethyl acetate (ester standard), and product in four separate test tubes. Add one drop of the sample to be tested to 1 mL of 0.5 M hydroxylamine hydrochloride (NH₂OH-HCl) in 95% ethanol in a test tube. Add 0.2 mL of a 6 M NaOH solution drop-wise and a boiling chip. Bring the mixture to a boil by heating in a water bath. Let the system cool and add 2 M HCl drop-wise until the pH is 2-3. If cloudiness develops, add 2 mL of 95% ethanol. Add 2 drops of 3% ferric chloride solution. A red-violet color is a positive test.

Analysis: IR and NMR spectroscopy spectroscopy

"The [crude or purified] product was analyzed by IR."

Analyze the IR spectrum of the alcohol (provided in lab, posted on Canvas). Obtain the IR spectrum of your product using NaCl plates and identify ester peaks. *Is there an OH peak in the product?*

Interpret the ¹H NMR spectra of <u>both ester products</u> (may be provided in lab, also posted on Canvas). Assign every hydrogen on the structure to a signal on the spectrum (integration, splitting, and expected & observed chemical shift).

Table	1. Clean-up & Safety	
-		

Clean-up	Safety
Liquid waste: aqueous layers and solutions from	H ₂ SO ₄ , HCl, hydroxamic acid, NaOH, and acetic acid
chemical test	are corrosive
Solid waste: Na ₂ SO ₄ , pipets, filter pipets	Acetic acid, ethyl acetate, and ethanol are flammable
After analysis, dispose of your product in the liquid	Isoamyl alcohol is an <i>irritant</i>
waste using a very small amount of ethanol from a	
wash bottle to aid the transfer.	
Wash all glassware and wipe down counters; return	Wear gloves & goggles throughout the experiment.
shared glassware to reagent counter	
Clean IR plates with acetone saturated with NaCl.	Be careful not to burn or melt the water hoses on the
Return plates to the desiccator after use.	hotplates!
Clean IR plates with acetone saturated with NaCl. Return plates to the desiccator after use.	Be careful not to burn or melt the water hoses on the hotplates!

All used pipets & broken glass go in the glass waste box. Please do not throw away glass in the trash as it creates an unexpected occupational hazard for our custodial staff. Thank you for participating in community set up & clean up tasks ©

References and Supplemental Reading

Palleros, D. R. "Preparation of Fruity Fragrances," *Experimental Organic Chemistry*, **2000**. Wiley: Hoboken.

Mohrig 4th edition: Chapter 7.1 (Reflux), 22.7-11 (NMR)

Klein 2nd edition: Chapter 15.1-6 (NMR), 20.10 & 20.15 (Fischer esterification)

McMurry 8th edition: Chapter 13.11 (¹H NMR splitting), 21.3 (Fischer Esterification), & 21.10 (¹H NMR of esters)

Follow Canvas Exp 4 Module...

Before Lab

- Read this PDF background, procedure, safety, pre-lab and in-lab questions
 - Option: listen to Caitlin read this document in the 8M Exp 4 Podcast
- Attend and/or watch lab lecture we go over everything you need for your assignments!
- Practice the lab online via Slugs@home platform sites.google.com/ucsc.edu/slugshome/home
- Complete the pre-lab questions at the end of this doc incorporated into Canvas quiz ©
 - Quiz due before your enrolled section check Canvas for due date
- Download the Exp 4 worksheet and prepare your lab notebook...

Lab Notebook Preparation – worksheet = template / outline to copy by hand into lab notebook

- Purpose: one-sentence summary of the main lab goals plus the reaction schemes
- Reagent Table add chemical properties; Wikipedia is a reliable source for chemical properties!
- Procedure with Diagrams complete before starting lab; sample on Canvas
 - Use the procedure on the previous pages to create your hand-drawn experimental instructions
 - Simple sketches & labels for all **equipment, chemical names** with **amounts**, & **transfers**
 - Format: Break it up with flow charts, bullet-points, comic strip, and/or whatever works for you!
 - Avoid copying the procedure word-for-word.
 - Make it easy for anyone to follow your procedure without referring to this document.
 - Slugs@home Exp 4 website Equipment & Safety pages; pictures & videos of the whole lab
 - The class notes include useful diagrams as well

During Lab

- Check the safety rules to dress for lab and arrive a few minutes early to Thimann Labs
- Pre-lab talk: tips for success and open Q&A
- Show your lab notebook pages to your TA
- Perform the experiment with a partner, fill out data & observations in **lab notebook**

After Lab – each partner submits separate, individual assignments

- Upload <u>Notebook Pages</u> to Canvas by midnight on lab day graded on completeness / participation
- Complete & upload the Lab Report on GradeScope (GS) due date on Canvas
 - o In-lab questions & experimental methods see last page of this document

Pre-lab Questions / Quiz - see your class notes!

1. Why is the reaction mixture extracted with **sodium bicarbonate (NaHCO₃)** and **sodium chloride (NaCl)** solution? What role does each salt play?

2. Calculate the **mass (mg)** and **volume (mL)** of alcohols and acetic acid that will be mixed from the mmol given in the procedure. Include these values in the reagent table in your notebook.

3. Determine the **limiting reagent** in the reactions. Calculate the **theoretical yield** of both syntheses in mg. Recall that catalysts cannot be the limiting reagent.

4. How is Le Chatelier's Principle on equilibrium used to increase the success of the esterification reaction?

5. Based on the techniques you have learned thus far in the organic chemistry lab, what are two methods that could be used to **separate unreacted alcohol from the ester**? Briefly explain **why** each would be expected to work. *Hint: you learned one of the techniques earlier this quarter; the other you learned in 8L.*

Take the Exp 4 pre-lab quiz before your enrolled section – check Canvas for due date

- The quiz incorporates the questions below the questions may be reworded.
- Be prepared with your responses to the pre-lab questions *before* starting the quiz.
- There is a 20-minute time limit on the quiz and you get two attempts.
 - Make sure you have enough time to complete the quiz you can't save and come back later.
 - o If you choose to re-take the quiz, your grade will be the highest of the two attempts.

Though we encourage collaboration in this class, this is an individual quiz.

• The responses should be a product of your original work so that you are assessed on *your* understanding of the material.

Sharing your quiz or your responses in any format (screenshots, email, CHEGG, social media, text, carrier pigeon, etc.) is in violation of the UCSC academic integrity policy.

LAB REPORT

Canvas Modules > Experiment 4 Report for submission details

Upload to GradeScope (GS) - due date on Canvas

- "Select Pages" to correlate your responses to the GS outline ☺
- o Option to submit with ONE partner one person uploads then "Add Group Member"

A. In-Lab Questions - see your lecture notes!

1. Show the full arrow-pushing **mechanism** for the assigned ester (apple or banana), including charged reaction intermediates to account for all bonds broken and formed. Draw the **full structures of the starting alcohol and final product** (define "R groups" as use the abbreviation only in the intermediates).

2. Draw the chemical reactions that **sodium bicarbonate** facilitates in the reaction workup. The equations should include the gas formed. *Note: sodium bicarbonate participates in two different (though similar) reactions.*

3. Report the **yield (mg)** and calculate the **percent yield** of the assigned synthesis. Discuss 2-3 suspected sources of **product loss** (exact parts of the procedure, such as transfers between containers, when you most likely lost product).

4. Interpret the **IR of the assigned alcohol and product**. Include the functional group, bond, expected and observed absorbances (wavenumbers, cm⁻¹). Briefly discuss whether the reaction went to completion (or not).

5. Report and interpret the **hydroxamic acid test** results. Draw the **chemical reaction** that occurred with your product – no abbreviations. What do the chemical test results suggest about the **success** of the fruity fragrances synthesis?

6. Interpret the ¹H NMR spectra provided in lecture for *both* (a) banana oil and (b) sour apple. Spectra provided in lecture notes and posted online. *Caitlin made supplemental videos on NMR interpretation of these esters - linked in the Canvas assignment for this report* ©

• Re-create separate typed tables for each ester in your report, including re-drawn structures with each set of H's labeled (A, B, C, etc.).

B. Experimental Methods

Use the **bold headings within the Exp 4 procedure** to get an idea of the level of detail to include in the experimental methods section (you may use those exact words!). You will need to **fill in your own data and descriptions** in place of "…" Simply report whether the "presence of an ester was confirmed by the **hydroxamic acid test**" (no procedural details). **IR** is the only form of characterization to report, as you are not directly analyzing your sample by NMR.

Writing guidelines and **sample experimental methods** are available on Canvas. Remember the sample experimental contains way more information than is pertinent to CHEM 8M students!